

PIP4K2B Antibody



PACO11254

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PIP4K2B

Uniprot

P78356

Synonyms:

phosphatidylinositol-5-phosphate 4-kinase, type II, beta; PIP4K2B; PI5P4KB; PIP5K2B; PIP5KIIB; PIP5KIIBeta ;

Immunogen:

Human PIP4K2B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A