

PIP5K1A Antibody



PACO11255

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF, IP

Recommended dilutions:

Protein Background:

Gene ID:

PIP5K1A

Uniprot

Q99755

Synonyms:

phosphatidylinositol-4-phosphate 5-kinase, type I, alpha;

Immunogen:

Human PIP5K1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A