

PIP5K1B Antibody



PACO11256

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PIP5K1B

Uniprot

O14986

Synonyms:

phosphatidylinositol-4-phosphate 5-kinase, type I, beta; PIP5K1B; MSS4; STM7 ;

Immunogen:

Human PIP5K1B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A