

PIKFYVE Antibody



PACO11257

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PIKFYVE

Uniprot

Q9Y2I7

Synonyms:

phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III; PIP5K3; CFD; FAB1; KIAA0981; MGC40423; PIKFYVE; PIP5K ;

Immunogen:

Human PIP5K3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A