

PITRM1 Antibody



PACO11263

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PITRM1

Uniprot

Q5JRX3

Synonyms:

pitrilysin metallopeptidase 1; PITRM1; KIAA1104; MGC138192; MGC141929; MP1; PreP; hMP1 ;

Immunogen:

Human PITRM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A