

# PKD2 Antibody



PACO11267

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PKD2

**Uniprot**

Q13563

**Synonyms:**

polycystic kidney disease 2 (autosomal dominant); PKD2; APKD2; MGC138466; MGC138468; PC2; PKD4 ;

**Immunogen:**

Human PKD2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A