

# PKIB Antibody



PACO11269

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PKIB

**Uniprot**

Q9C010

**Synonyms:**

protein kinase (cAMP-dependent, catalytic) inhibitor beta; PKIB; FLJ23817; PRKACN2 ;

**Immunogen:**

Human PKIB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A