

PKP4 Antibody



PACO11277

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PKP4

Uniprot

Q99569

Synonyms:

plakophilin 4; PKP4; FLJ31261; FLJ42243; p0071 ;

Immunogen:

Human PKP4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A