

# PLA2G4E Antibody



PACO11280

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PLA2G4E

**Uniprot**

Q3MJ16

**Synonyms:**

phospholipase A2, group IVE; PLA2G4E; FLJ45651; MGC126633; MGC126661 ;

**Immunogen:**

Human PLA2G4E.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A