

# PLCD4 Antibody



PACO11288

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PLCD4

**Uniprot**

Q9BRC7

**Synonyms:**

phospholipase C, delta 4; PLCD4; MGC12837 ;

**Immunogen:**

Human PLCD4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A