

PLEKHA5 Antibody



PACO11296

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PLEKHA5

Uniprot

Q9HAU0

Synonyms:

pleckstrin homology domain containing, family A member 5; PLEKHA5; FLJ10667; FLJ31492; KIAA1686; PEPP2 ;

Immunogen:

Human PLEKHA5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A