

# PLEKHA8 Antibody



PACO11297

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PLEKHA8

**Uniprot**

Q96JA3

**Synonyms:**

pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 8; PLEKHA8; FAPP2; MGC3358 ;

**Immunogen:**

Human PLEKHA8.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A