

# PLVAP Antibody



PACO11312

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PLVAP

**Uniprot**

Q9BX97

**Synonyms:**

plasmalemma vesicle associated protein; PLVAP; FELS; PV-1; PV1; gp68 ;

**Immunogen:**

Human PLVAP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A