

PMFBP1 Antibody



PACO11316

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PMFBP1

Uniprot

Q8TBY8

Synonyms:

polyamine modulated factor 1 binding protein 1; PMFBP1; DKFZp434G131; FLJ40146 ;

Immunogen:

Human PMFBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A