

PMM2 Antibody



PACO11318

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PMM2

Uniprot

O15305

Synonyms:

phosphomannomutase 2; PMM2; CDG1; CDG1a; CDGS ;

Immunogen:

Human PMM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A