

PMP2 Antibody



PACO11319

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PMP2

Uniprot

P02689

Synonyms:

peripheral myelin protein 2; PMP2; FABP8; M-FABP; MP2; P2 ;

Immunogen:

Human PMP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A