

PMPCB Antibody



PACO11320

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PMPCB

Uniprot

O75439

Synonyms:

peptidase (mitochondrial processing) beta; PMPCB; Beta-MPP; MPP11; MPPB; MPPP52 ;

Immunogen:

Human PMPCB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A