

# PNRC2 Antibody



PACO11335

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PNRC2

**Uniprot**

Q9NPJ4

**Synonyms:**

proline-rich nuclear receptor coactivator 2; PNRC2; FLJ20312; MGC99541 ;

**Immunogen:**

Human PNRC2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A