

# POF1B Antibody



PACO11337

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POF1B

**Uniprot**

Q8WVV4

**Synonyms:**

premature ovarian failure, 1B; POF1B; FLJ22792; POF ;

**Immunogen:**

Human POF1B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A