

POFUT2 Antibody



PACO11339

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

POFUT2

Uniprot

Q9Y2G5

Synonyms:

protein O-fucosyltransferase 2; POFUT2; C21orf80; FUT13 ;

Immunogen:

Human POFUT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A