

POLDIP2 Antibody



PACO11344

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

POLDIP2

Uniprot

Q9Y2S7

Synonyms:

polymerase (DNA-directed), delta interacting protein 2; POLDIP2; DKFZp586F1524; PDIP38; POLD4 ;

Immunogen:

Human POLDIP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A