

POLDIP3 Antibody



PACO11345

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

POLDIP3

Uniprot

Q9BY77

Synonyms:

polymerase (DNA-directed), delta interacting protein 3; POLDIP3; KIAA1649; PDIP46; SKAR ;

Immunogen:

Human POLDIP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A