

POLG2 Antibody



PACO11347

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

POLG2

Uniprot

Q9UHN1

Synonyms:

polymerase (DNA directed), gamma 2, accessory subunit; POLG2; HP55; MTPOLB; PEOA4; POLB; POLG-BETA; POLGB ;

Immunogen:

Human POLG2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A