

# POLR2L Antibody



PACO11361

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POLR2L

**Uniprot**

P62875

**Synonyms:**

polymerase (RNA) II (DNA directed) polypeptide L, 7.6kDa; POLR2L; RBP10; RPABC5; RPB10; RPB10beta; RPB7.6; hRPB7.6; hsRPB10b ;

**Immunogen:**

Human POLR2L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A