

# POLR3D Antibody



PACO11364

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POLR3D

**Uniprot**

P05423

**Synonyms:**

polymerase (RNA) III (DNA directed) polypeptide D, 44kDa; POLR3D; BN51T; RPC4; RPC53; TSBN51 ;

**Immunogen:**

Human POLR3D.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A