

PACO11366

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

POLR3K

Uniprot

Q9Y2Y1

Synonyms:

polymerase (RNA) III (DNA directed) polypeptide K, 12.3 kDa; POLR3K; C11; C11-RNP3; My010; RPC10; RPC11; hRPC11 ;

Immunogen:

Human POLR3K.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A