

# POMT2 Antibody



PACO11368

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POMT2

**Uniprot**

Q9UKY4

**Synonyms:**

protein-O-mannosyltransferase 2; POMT2; DKFZp686G10254; FLJ22309 ;

**Immunogen:**

Human POMT2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A