

POPDC3 Antibody



PACO11374

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

POPDC3

Uniprot

Q9HBV1

Synonyms:

popeye domain containing 3; POPDC3; MGC22671; POP3; RP11-99L11.2; bA355M14.1 ;

Immunogen:

Human POPDC3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A