

PPEF1 Antibody



PACO11391

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPEF1

Uniprot

O14829

Synonyms:

protein phosphatase, EF-hand calcium binding domain 1; PPEF1; PP7; PPEF; PPP7C ;

Immunogen:

Human PPEF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A