

# PPIL3 Antibody



PACO11403

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PPIL3

**Uniprot**

Q9H2H8

**Synonyms:**

peptidylprolyl isomerase (cyclophilin)-like 3; PPIL3; CYPJ ;

**Immunogen:**

Human PPIL3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A