

PPIL6 Antibody



PACO11406

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPIL6

Uniprot

Q8IXY8

Synonyms:

peptidylprolyl isomerase (cyclophilin)-like 6; PPIL6; MGC27054; MGC27056; MGC41939; bA425D10.6; dJ919F19.1 ;

Immunogen:

Human PPIL6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A