

PPM1B Antibody



PACO11408

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPM1B

Uniprot

O75688

Synonyms:

protein phosphatase 1B (formerly 2C), magnesium-dependent, beta isoform; PPM1B; MGC21657; PP2C-beta-X; PP2CB; PP2CBETA; PPC2BETAX ;

Immunogen:

Human PPM1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A