

PPP1R2P9 Antibody



PACO11424

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPP1R2P9

Uniprot

O14990

Synonyms:

protein phosphatase 1, regulatory (inhibitor) subunit 2 pseudogene 9; PPP1R2P9; I-4 ;

Immunogen:

Human PPP1R2P9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A