

PPP1R9B Antibody



PACO11428

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPP1R9B

Uniprot

Q96SB3

Synonyms:

protein phosphatase 1, regulatory (inhibitor) subunit 9B; PPP1R9B; FLJ30345; PPP1R6; PPP1R9; SPINO; Spn ;

Immunogen:

Human PPP1R9B.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A