

PPP2R2C Antibody



PACO11434

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PPP2R2C

Uniprot

Q9Y2T4

Synonyms:

protein phosphatase 2 (formerly 2A), regulatory subunit B, gamma isoform; PPP2R2C; B55-GAMMA; IMYPNO; IMYPNO1; MGC33570; PR52; PR55G ;

Immunogen:

Human PPP2R2C.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A