

# PPP3CC Antibody



PACO11437

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PPP3CC

**Uniprot**

P48454

**Synonyms:**

protein phosphatase 3 (formerly 2B), catalytic subunit, gamma isoform; PPP3CC; CALNA3 ;

**Immunogen:**

Human PPP3CC.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A