

# PRDM4 Antibody



PACO11450

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRDM4

**Uniprot**

Q9UKN5

**Synonyms:**

PR domain containing 4; PRDM4; MGC45046; PFM1 ;

**Immunogen:**

Human PRDM4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A