

PRELID1 Antibody



PACO11456

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PRELID1

Uniprot

Q9Y255

Synonyms:

PRELI domain containing 1; PRELID1; CGI-106; MGC87972; PRELI; PX19 ;

Immunogen:

Human PRELID1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A