

# PRKAA2 Antibody



PACO11464

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRKAA2

**Uniprot**

P54646

**Synonyms:**

protein kinase, AMP-activated, alpha 2 catalytic subunit; PRKAA2; AMPK; AMPK2; PRKAA ;

**Immunogen:**

Human PRKAA2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A