

PRPF40B Antibody



PACO11499

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PRPF40B

Uniprot

Q6NWX9

Synonyms:

PRP40 pre-mRNA processing factor 40 homolog B (S. cerevisiae); PRPF40B; HYPC ;

Immunogen:

Human PRPF40B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A