

PRPS1L1 Antibody



PACO11504

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PRPS1L1

Uniprot

P21108

Synonyms:

phosphoribosyl pyrophosphate synthetase 1-like 1; PRPS1L1; PRPS1; PRPS3; PRPSL; PRS-III ;

Immunogen:

Human PRPS1L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A