

# PRTFDC1 Antibody



PACO11511

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRTFDC1

**Uniprot**

Q9NRG1

**Synonyms:**

phosphoribosyl transferase domain containing 1; PRTFDC1; FLJ11888; HHGP ;

**Immunogen:**

Human PRTFDC1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A