

PSAPL1 Antibody



PACO11515

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PSAPL1

Uniprot

Q6NUJ1

Synonyms:

prosaposin-like 1; PSAPL1; FLJ40379 ;

Immunogen:

Human PSAPL1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A