

PSMC3IP Antibody



PACO11542

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PSMC3IP

Uniprot

Q9P2W1

Synonyms:

PSMC3 interacting protein; PSMC3IP; GT198; HOP2; HUMGT198A; TBPIP ;

Immunogen:

Human PSMC3IP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A