

# PSMD11 Antibody



PACO11546

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PSMD11

**Uniprot**

O00231

**Synonyms:**

proteasome (prosome, macropain) 26S subunit, non-ATPase, 11; PSMD11; MGC3844; Rpn6; S9; p44.5 ;

**Immunogen:**

Human PSMD11.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A