

PTPRM Antibody



PACO11598

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PTPRM

Uniprot

P28827

Synonyms:

protein tyrosine phosphatase, receptor type, M; PTPRM; MGC166994; PTPRL1; R-PTP-MU; RPTPM; RPTPU; hR-PTPu ;

Immunogen:

Human PTPRM.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A