

RAB11FIP2 Antibody



PACO11631

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

RAB11FIP2

Uniprot

Q7L804

Synonyms:

RAB11 family interacting protein 2 (class I); RAB11FIP2; KIAA0941; Rab11-FIP2; nRip11 ;

Immunogen:

Human RAB11FIP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A