

RAB37 Antibody



PACO11652

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RAB37

Uniprot

Q96AX2

Synonyms:

RAB37, member RAS oncogene family; RAB37; FLJ30284; FLJ32507 ;

Immunogen:

Human RAB37.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A