

RAB39B Antibody



PACO11654

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RAB39B

Uniprot

Q96DA2

Synonyms:

RAB39B, member RAS oncogene family;

Immunogen:

Human RAB39B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A