

RAB3IL1 Antibody



PACO11660

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

RAB3IL1

Uniprot

Q8TBN0

Synonyms:

RAB3A interacting protein (rabin3)-like 1; RAB3IL1; GRAB ;

Immunogen:

Human RAB3IL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A